

LIVESTOCK WATER SUITABILITY



CLIENT INFORMATION

Client Name _____ Account # _____

Submitted For _____

Address _____

City _____ State _____ ZIP _____

Email _____ Phone # _____

SAMPLE INFORMATION

Sample ID _____ Collection Date ____/____/____ Collection Time _____

Livestock Suitability:	Livestock Plus:	Optional Additional Tests: <i>(requires additional bottle)</i>
pH, Conductivity, Total Dissolved Solids, Nitrate, Sulfate, Chloride, Sodium, Iron	Livestock Suitability + Calcium, Magnesium, Hardness, Alkalinity, Alkalinity Bicarbonate, S.A.R., Manganese	Total/E.coli Coliform Bacteria

SUBMITTING A SAMPLE

- 1 Use the sampling containers provided by AgSource Laboratories. ***Some sample containers contain a liquid acid or powdered preservative. Use caution when collecting sample and do not dispose of preservative.**
- 2 Carefully open the bottle and take care not to contaminate the inside of the bottle or lid.
- 3 Place bottle under a slowly flowing stream of water and fill completely. It is important to supply an adequate amount of sample to get a representative measurement. This also ensures that the laboratory has enough sample to run the standard tests and Quality Assurance samples.
- 4 Clearly fill out all sample tags and submittal form *completely*.
- 5 Ship your samples directly to:
AgSource Laboratories
1701 Detroit Street, PO Box 247
Ellsworth, IA 50075

***Samples must arrive at the laboratory within 48 hours of collection. Samples accepted Monday - Thursday.**